


HySafe Qui-Connect


Modular Forceps for GYN Hysteroscopy

360° rotatable, flushable, easy click, rigid

Insert with sheath	Insert with sheath			Fr.	WL (mm)	Handle with overload protection	Complete instrument
	Insert with sheath	Fr.	WL (mm)				
Grasping forceps with alligator jaws	8642.6508	5	340	8083.001	Types	<ul style="list-style-type: none"> • Insert with sheath • Handle 	
	8992.6308	7	310				
Biopsy forceps, -spoon	8642.6008	5	340				
Biopsy punch	8642.6318	5	340				
Micro scissors	8642.6818	5	340				
Hooked scissors	8992.6408	7	310				

Technical specifications subject to change without notice.

639-01.01-0408USA



HySafe Qui-Connect

Modular Forceps for GYN Hysteroscopy

- Easy, quick assembly
- Modular forceps and scissors
- Flushable
- Spring overload protection




For more information, contact your local Richard Wolf representative or call Customer Service at (800) 323-WOLF (9653).

QUI-CONNECT MODULAR FORCEPS

THE ADVANTAGES – AT A GLANCE!

High efficiency:



Economical due to modularity

High efficiency:

- Due to overload protection, the service life of the jaws is significantly longer
- All the jaws can be implemented using just one forceps handle

Hygiene:

- Axial flushing possible for optimal cleaning

Durability:

- Optimum durability due to overload protection

Excellent handling:

- User-friendly disassembly using the Qui-Connect mechanism
- No adverse effect due to integrated flush connection
- User friendly due to marking aid
- Jaws rotatable through 360°
- Simple assembly in any position
- Jaws with easy, quick, replacement

For many years, Richard Wolf has ranked amongst the most successful innovators throughout the world who have been actively engaged in the research and development of medical products and product systems. Since the demands for efficient application solutions are continually on the increase, Richard Wolf is constantly endeavoring to offer new, unique products for applications that are already available, or for new ones.

By constructing this modular design with its new Qui-Connect mechanism, in combination with optimized durable components, its concept has been implemented in an excellent way. Due to the new system's overload protection a service life is guaranteed of at least three times longer.

With this new development, important details in regards to efficiency, user-friendliness and hygiene have been greatly improved. Interdisciplinary implementation is possible in urology, gynecology and surgery.

Parallel to the rapid development and implementation of new application possibilities we have addressed our customers' wish for developing the new Qui-Connect Forceps with flushable forceps and scissors.

Fast easy assembly
Thanks to HySafe Qui-Connect; as simple as turning over this page. See and feel the difference!

Hygiene:



Thorough cleaning of the jaw



and interior lumen

Durability:



Spring overload protection



provides a significantly longer service life

Excellent handling:



Simplest assembly and disassembly



using the Qui-Connect mechanism

